FISĂ TEHNICĂ 34077

ROMANA

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Se aplica garantia standard Gima B2B de 12 luni.

GIMA 34077

Zhangjiagang Tengda Machinery Manufacturer Co.,Ltd
No.999 Hongqi Road, Changyinsha Modern Agricultural Demonstration Area,
Zhangjiagang City, Jiangsu, 215623, China Made in China
SUNGO Europe B.V., Fascinatio Boulevard 522, Unit 1.7, 2909VA
Capelle aan den Ijssel, The Netherlands
Gima S.p.A., Via Marconi, 1 - 20060 Gessate (MI) Italy

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

II. Technical Features

1. Product size (L x W x H): 188 x 53 x 21 cm
2. Folded size (L x W x H): 92 x 50 x 11 cm
3. Packing size (2pcs/carton): 96 x 54 x 21 cm
4. Net Weight: 8 kg Gross Weight: 9 kg
5. Load bearing: less than 159 kg

II. How to use

This separable first-aid folding stretcher is mainly used for transporting patients with fractures. Both ends of the stretcher have hinge-type clutches to separate the stretcher into a right and left part. This allows for the patient to be quickly and gently placed onto the stretcher, operation table or bed.

III. Applied for

Approved by clinic practice, this item is suitable for ambulance, hospital battlefield gym, to carry the wounded & patients. It also can be used in general operation.

IV. Cautions

1. When carrying patients with this stretcher, be sure lock tightly the fender locking device & fasten the safety belt for safety sake.
2. In general operation, be careful not to out torn the stretcher surface.

V. Maintenance

1. Keep clean as routine (incl. sterilization).
2. Often check if any parts loose or not.

VI. Storage & transportation

1. Store this product in the place of damp proof & non-corrosion.
2. Common transporting vehicle can carry this product.

VII. Package

1. Packed in high durable carton with wooden base.
2. Packing: certificate of inspection & user's manual & packing list.

	Manufacturer		Keep in a cool, dry place	CE	Medical Device compliant with Regulation (EU) 2017/745
REF	Product code		Keep away from sunlight		Caution: read instructions (warnings) carefully
LOT	Lot number	MD	Medical Device	EC REP	Authorized representative in the European community
	Consult instructions for use		Date of manufacture	UDI	Unique device identifier
	Imported by				

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

GIMA 34077

Zhangjiagang Tengda Machinery Manufacturer Co.,Ltd
No.999 Hongqi Road, Changyinsha
Modern Agricultural Demonstration Area,
Zhangjiagang City, Jiangsu, 215623, China
Made in China

REF **TD010112**

EC REP SUNGO Europe B.V.
Fascinatio Boulevard 522, Unit 1.7, 2909VA
Capelle aan den Ijssel, The Netherlands

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

134077-M-Rev.7-04.25