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RO

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EN

Instructions for use
RAISED TOILET SEAT WITH LOCK

REF

NORMATIVE

This product is a Class I medical device and complies with Regulation (EU) 2017/745 of the European Parliament and of the Council of April 5, 2017 on Medical Devices.

This product should not be disposed of in the waste container. There are specific collection points where you can deposit them. For the disposal of the container and the product, strictly comply with the legal regulations of your community.

INTENDED USE

Support product for lifting fixed to the toilet. It provides stability and facilitates the evacuation functions of people with disabilities or reduced mobility. Its use is also advised after hip surgery.

AVAILABLE SIZES

WITHT LID

- 10 cm (4 inch)
- 15 cm (6 inch)

TECHNICAL CHARACTERISTICS

- Material: rigid plastic
- Measures: 36 (width) x 40 cm (depth)
- Height: depends on the chosen model
- Net Weight depeding on the model: 0,81 – 1,44 Kg
- Max. Weight capacity: 120 Kg

PRECAUTIONS

- Before each use, you must make sure that the product is in perfect conditions.
- Do not use the product in case of a visible defect. In case of missing, ill-fitting or damaged parts, ask your pharmacist or distributor for advice.
- Ensure that the toilet lift is fixed with the screws in the toilet bowl.
- If instructions, warnings, and notes in the user manual are not followed, it may cause physical injury to you or others, or damage to the product.
- Do not make any modifications (mechanical alterations) to the item so as not to compromise its safety of use.
- Do not exceed the maximum weight.

CONTRAINDICATIONS

The use of the raised toilet is not recommended for people with strong perception disorders, balance disorders or inability to sit.

CARE & MAINTENANCE RECOMMENDATIONS

- Clean and dry periodically. Remove dirt with a soft, damp cloth and neutral soap. Dry later. Do not use aggressive or abrasive chemical cleaners; they can damage the material of the product. Also, do not use greasy products that can degrade materials or make it slippery.
- When not in use, store the product in a dry place at room temperature.

LIFE SPAN

The expected life span of this product is 2 years under normal conditions of use, safety and maintenance.

Beyond this period, the product can be used as long as it is in perfect conditions of use and maintenance. In any case, the product may require the replacement of certain parts or components due to normal wear in a period shorter than indicated. The useful life is influenced by its conditions of use and storage and by the regularity of its maintenance.

WARRANTY

The warranty is for a period of three years (6 months on batteries) from the date of sales invoice, as long as the product has not been altered in its original configuration. Within the warranty period we will remedy any manufacturing defect free of charge, either by repairing, replacing parts or providing a new product, at our discretion. To opt for repair or replacement, the owner of the product must go to the retail establishment, which will be the only one in charge of managing the warranty with the manufacturer.

This warranty will not be valid in cases of improper use or handling of the product, deterioration of the product due to external agents such as water or other harmful chemicals, obstructive or corrosive products, due to inadequate maintenance, lack of cleaning, component loss or the use of non-original spare parts of our brand. For more information check our general conditions of sale.

WARNING

- If you observe any deficiency, anomaly or failure, immediately notify the vending establishment or distributor to proceed with its change or repair.
- Any serious incident related to the product must be reported to TESA MED SOLUTIONS and to the competent authority of the Member State in which the user and/or the patient are established.